

Cyclone Collector



- 1 From model trimmer or sandblaster
- 2 Cyclone collector
- 3 Debris pail
- 4 Outlet to dust collector
- 5 Dust collector

SANDBLASTING · DIE-STONE · ZIRCONIA · DRY MODEL TRIMMING · ARCH TRIMMING

AA249 fits 1.5" or 2" hose and includes:

DIRTY AIR INLET FROM SOURCE

- (1) rubber fitting for 2" diameter hose (#H481)
 - (1) rubber fitting for 1.5" diameter hose (#H842)
- If using 2" (50mm) diameter hose, connect #H481 to hose that you are using
 If using 1.5" (38mm) diameter hose, connect #H842 to hose that you are using

CLEAN AIR OUTLET TO DUST COLLECTOR

- 2" diameter hose, 6 ft long (#H154-06)
 - (2) rubber fittings (#H481)
 - (1) rubber fitting (#H767)
- If dust collector includes 2" (50mm) diameter inlet, connect two (2) #H481 fittings to hose supplied
 If dust collector includes 1.5" (38mm) diameter inlet, connect #H481 and #H767 fittings to hose supplied



The Cyclone Collector separates the larger particles from the air stream increasing the efficiency of the filtration system by preventing frequent clogging and loss of suction.

Cyclone Collector's small footprint allows it to fit under any workbench making it the perfect choice for dental labs.

Manufactured using static dissipating polymers, the Cyclone Collector offers additional static reduction for dental applications.

Part #: AA249
 Dimensions: 27" h (686mm)* x 12" dia (305mm)
 *add 4" (102mm) for hose
 Capacity: 0.8 ft³ (5 gal / 23 l)